

Gas Chromatograph Mass Spectrometer

GCMS-TQ8040

UFMS
ULTRA FAST MASS SPECTROMETRY





GCMS-TQ8040 *with Smart Technologies*

Smart Performance, Smart Productivity, Smart Operation

Finally, a triple quadrupole GCMS *Smart* enough for everyday use in your laboratory.




The Shimadzu GCMS-TQ8040 is the first triple quadrupole with **Smart Performance** for low detection limits and Scan/MRM, **Smart Productivity** for high efficiency sample throughput, and **Smart Operation** for quick and easy method development. These three smart technologies contribute to **Smart MRM**, and provide the most accurate, cost-effective, and easy-to-use triple quadrupole GCMS you have ever imagined.

Smart Performance

The ion source, incorporating a highly efficient collision cell, enables high-sensitivity and long-term stable analyses. In addition, with high-speed scan control technology called Advanced Scanning Speed Protocol (ASSP), multiple analysis modes can be measured simultaneously, enabling a great deal of information to be obtained in a single analysis.

> | P. 6 to **High Sensitivity and Long-Term Stability**



Smart Productivity

Simultaneous multi-component analysis is now possible for hundreds of target compounds, dramatically improving productivity. Smart MRM produces GC-MS/MS methods with up to 32,768 transitions in a single run.

> | P. 8 to **Smart Productivity**



Smart Performance Productivity Operation

Smart Operation

GC-MS/MS analysis requires multiple settings that can be confusing to the average operator. With **Smart MRM**, the GCMS-TQ8040 software sets the analytical conditions automatically, making method development painless, fast, and easy.

> | P. 10 to **Smart Operation**





Smart Performance

High Sensitivity and Long-Term Stability

■ Ion Source Featuring High Sensitivity and Long-Term Stability



The filament and the ion source box have been separated, which reduces the impact of the filament potential on the interior of the ion source. In addition, a shield is provided to block radiant heat generated by the filament, achieving a uniform temperature inside

the ion source box. Thus, active spots inside the ion source are not prone to occur, enabling high-sensitivity analysis with long-term stability. (Patent: US7939810)




■ High-Performance Quadrupole Mass Filter

A high-accuracy mass filter equipped with a pre-rod, as well as patented electric field control technology, result in high-accuracy mass separation performance. The pre-rods also minimize contamination of the quadrupole, eliminating the need for quadrupole maintenance.




Response of a Pesticide (Fenitrothion) in a Sesame Extract



■ OFF-AXIS Ion Optics


Lower detection limits are achieved by OFF-AXIS Ion Optics (Patent Pending). Meta-stable and neutral ions are removed without sacrificing sensitivity. Helium buffer gas is not required in the CID cell.



■ High-Efficiency Collision Cell UFSweeper

Shimadzu's proprietary UFSweeper technology achieves high-speed MRM analysis at speeds up to 800 transitions per second. It sweeps residual ions from the collision cell to provide high-efficiency CID and fast ion transport. Rapid ion removal minimizes cross-talk and enables trace analysis (patent pending).

UFMS
ULTRA FAST MASS SPECTROMETRY




Analysis of Residual Pesticides (Isoprothiolane 1 pg/μL)

■ Sensitivity and Repeatability in Single GC/MS Mode


The high-efficiency ion source provides the foundation of an ion generation and transmission system, which creates and then delivers ions to the detector, resulting in a GC/MS with the maximum possible sensitivity and repeatability. These features

are not realized just for MRM measurements by GC-MS/MS, but also for scan and SIM measurements in single quadrupole modes, even with the most reactive compounds.


Q1 with post-rod achieves outstanding ion transfer efficiency




UFSweeper provides high ion transfer efficiency



Thiobencarb 5 pg/μL
Left: GCMS-QP2010 Ultra, Right: GCMS-TQ8040
(for RSD determination, n = 5)



Scan Mass Spectrum of Fenitrothion



Mass Spectrum Registered in NIST Library


Smart Productivity

▮ Innovative Technologies that Improve Accuracy and Throughput


UFMS in the Multiple Reaction Mode can acquire over 800 transitions per second and over 32,000 transitions in a single analysis. The **Smart MRM** technology automatically adjusts the analytical dwell time for each transition, only acquiring data during peak elution, to fully optimize sensitivity. For example, analysis of more than 400 pesticides that used to require two or three methods can now be accomplished in a single acquisition method created by **Smart MRM**, significantly increasing laboratory throughput.

The UFMS technology guarantees a minimum of ten data points across each peak for optimum sensitivity and repeatability. Highly accurate, low-level detection of multiple components in complex matrices is finally possible by UFMS and **Smart MRM**.


> | P. 11 to **Smart MRM**



Simultaneous Analysis of 439 Pesticides Using UFMS and **Smart MRM**




Mass Chromatogram and %RSD (5 pg/μL)



Reduce Analysis Costs and Minimize Downtime


If simultaneous, multi-component analysis can be performed in 1/2 to 1/3 the time of existing systems, then two to three times the number of samples can be analyzed in the same period of time, and return on investment is improved. In addition, the frequency of maintenance, such as replacing glass liners and columns, is reduced, thereby minimizing downtime.



Twin Line MS System Eliminates the Need to Vent the MS

The GCMS-TQ8040 is capable of accepting installation of two narrow-bore capillary columns into the MS simultaneously. This allows you to switch applications without venting the MS.

Simply decide which column is best for your analysis and choose the associated injection port.



Methamidophos in Ginger (10 pg/μL),
Analyzed on Two Dissimilar Columns Using **Smart MRM**

CID gas control is a method parameter, allowing acquisition of GC-MS and GC-MS/MS data in the same batch. By coupling this with the Twin Line MS System, analysis of phenols by SIM and analysis of pesticides by MRM are possible in a single batch without venting the MS.

	Method File	Data File	
1	Phenols_SIM.qgm	DATA_01.qgd	SIM
2	Phenols_SIM.qgm	DATA_02.qgd	
3	Phenols_SIM.qgm	DATA_03.qgd	
4	Pesticide_MRM.qgm	DATA_04.qgd	MRM
5	Pesticide_MRM.qgm	DATA_05.qgd	

SIM
GC-MS (CID gas—Off)

MRM
GC-MS/MS (CID gas—On)

Run SQ and TQ Methods in a Single Sequence

*) The Twin Line MS System requires an optional installation kit and is limited to certain column dimensions.

Smart Operation

▮ Simplified GC-MS/MS method development

In order to create an MRM analysis method, complex parameter settings are required, such as for the optimum transitions and collision energies.

Smart MRM makes method development quick and easy.

Whether starting from scratch to optimize transitions and collision energies for new compounds, or starting from an MRM database of known target analytes to build a custom MRM method,


Smart MRM takes the stress and difficulty out of method development.

MRM Optimization Tool

■ Optimize MRM Transitions Automatically

Determining and optimizing MRM transitions for new compounds can require significant development time. The “MRM Optimization Tool” automates the process by collecting product ion scan data and finding the optimum collision energy for each transition.

Once established, the transitions are registered to one of the Shimadzu “Smart Database” files, and the MRM or Scan/MRM methods are created using **Smart MRM**.



Smart Database

■ Compounds and Optimized Transitions

The Shimadzu “Smart Database” is a database of related compounds (e.g. pesticides, drugs, metabolites, etc.) with optimized transitions and collision energies, CAS registry numbers, and Retention Indices (RI). The user can select from hundreds of pre-registered compounds in one of the “Smart Database” files, or add their own optimized transitions. The user selects the compounds to be analyzed, and

Smart MRM builds the MRM or Scan/MRM acquisition method from the “Smart Database” with the push of a button.

In addition to MRM information, the database file can contain Scan and SIM ion information, mass spectra and calibration curve information from the internal standard method. This allows users to easily create their own database.

Serial#	Type	Acq. Mode	ISTD Group	Level1 Conc (IS)	Method No.	Compound Name (E)	Ret. Index 1	Cas#	Ion1			
							Method1		Type	m/z	CE	Rati
1	Target	MRM			1	Aldicarb deg.	881	0 - 00 - 0	T	115.1>68.0	8	100.00
2	Target	MRM			1	DCIP	1058	108 - 60 - 1	T	121.1>45.0	4	100.00
3	Target	MRM			1	Aldoxycarb deg.	1135	0 - 00 - 0	T	80.0>65.0	6	100.00
4	Target	MRM			1	Chlofentezine deg.	1182	0 - 00 - 0	T	137.0>102.0	14	100.00
5	Target	MRM			1	Hymexazol	1196	10004 - 44 - 1	T	99.0>71.0	8	100.00
6	Target	MRM			1	Methamidophos	1236	10265 - 92 - 6	T	141.0>95.0	8	100.00
7	Target	MRM			1	Dichlorvos	1263	62 - 73 - 7	T	109.0>79.0	8	100.00
8	Target	MRM			1	Nereistoxin	1283	0 - 00 - 0	T	149.1>71.1	8	100.00
9	Target	MRM			1	Allidochlor	1296	93 - 71 - 0	T	132.1>56.0	8	100.00
10	Target	MRM			1	Dichlobenil	1358	1194 - 65 - 6	T	170.9>136.0	14	100.00
11	Target	MRM			1	EPTC	1364	759 - 94 - 4	T	189.1>128.1	4	100.00
12	Target	MRM			1	Blphenyl	1394	92 - 52 - 4	T	154.1>128.1	22	100.00
13	Target	MRM			1	Propamocarb	1398	24579 - 73 - 5	T	188.2>72.0	4	100.00

Smart MRM



■ Automatic Method Creation

The Smart MRM technology automatically creates methods with measurement times optimized for each component based on the Smart database. The Automatic Adjustment of Retention Time (AART) function incorporated in the system estimates retention times with high accuracy. When creating methods for simultaneous multicomponent analysis, the complicated process of configuring

measurement parameters made it difficult to prepare appropriate methods. By using the Smart MRM function, however, it is possible to automatically create methods in which data are acquired with high sensitivity only during the elution time of the target components. In addition to MRM methods, SIM methods can be created.

The screenshot shows the 'Create Method File' dialog box with 'Instrument Type' set to 'TO Series' and 'Lang' set to 'Lamp'. Below it, the 'Parameter' section shows 'Ret. Index for AART' set to 'Ret. Index 1' and 'Template Method File' set to 'C:\GCMS\SmartMRM\SmartMRM27'. The 'MS Table Parameter' dialog box is also visible, showing 'MRM, SIM Parameter' with 'Loop Time (MRM, SIM)' set to '0.30 sec' and 'Required Processing Time: R.T.s' set to '0.30 min'. The 'Scan Mode' is set to 'ON' and 'Scan Parameter' is set to 'Event Time of Scan: 0.10 sec'.

Parameter Setting


■ Automatic Adjustment of Compound Retention Time (AART) (Automatic Adjustment of Retention Time)

Compound Name (E)	Ret. Time
Permethrin ethyl	23.286
Isopropalin	23.415
Isodrin	23.436
Cyprodinil	23.458
Isopropalin-methyl	23.500
MGK 264-2	23.565
(E)-Chlorfenvinphos	23.629
Dimethametryn	23.715
Metazachlor	23.737
Pendimethalin	23.780
Ethionazole	23.801
Penconazole	23.844
Heptachlor-epoxy-epoxide	23.887
oxy-Chlordane	23.930
Chlorzoxatone	23.952
(Z)-Pyrethrin	23.982
Thiabendazole	23.973
Allethrin-1,2	23.995
Phosfolan	24.016

Automatic Estimation of Retention Times

200 : Permethrin ethyl (3.3 ms)
201 : Isopropalin (3.3 ms)
202 : Isodrin (3.3 ms)
203 : Cyprodinil (3.3 ms)
204 : Isopropalin-methyl (3.5 ms)
205 : MGK 264-2 (3.3 ms)
206 : (E)-Chlorfenvinphos (3.3 ms)
207 : Dimethametryn (2.7 ms)
208 : Metazachlor (2.7 ms)
209 : Pendimethalin (2.7 ms)
210 : Ethionazole (2.7 ms)
211 : Aniazine (2.5 ms)
212 : Penconazole (2.7 ms)
213 : Heptachlor-epoxy-epoxide (2.7 ms)
214 : oxy-Chlordane (2.7 ms)
215 : Chlorzoxatone (2.7 ms)
216 : (Z)-Pyrethrin (2.7 ms)
217 : Thiabendazole (2.7 ms)
218 : Allethrin-1,2 (2.7 ms)
219 : Phosfolan (2.5 ms)
220 : Tolyfluand (2.7 ms)
221 : Heptachlor-endo-epoxide (2.7 ms)
222 : Bromfenoxim-methyl (2.7 ms)

Dwell Time Optimization



The AART function adjusts the retention times of target components based on linear retention indices (LRI) and the retention times of n-alkanes. The AART function easily adjusts acquisition and processing method parameters simultaneously.



Automatic Creation of Analysis Methods


Compound Name	Acq. Mode	Event Time(sec)	Ch1 m/z	Ch1 CE	Ch2 m/z	Ch2 CE	Ch3 m/z	Ch3 CE
44-14 MGK 264-2	MRM	0.018	164.10>98.00	12.00	111.10>82.00	8.00	164.10>67.00	8.00
44-15 (E)-Chlorfenvinphos	MRM	0.018	323.00>267.00	16.00	267.00>159.00	18.00	267.00>203.00	12.00
44-16 Dimethametryn	MRM	0.018	212.10>122.10	12.00	212.10>94.00	22.00	212.10>71.00	18.00
44-17 Metazachlor	MRM	0.018	209.10>132.10	18.00	133.10>117.10	24.00	211.10>132.10	20.00
45-1 Diphenamid	MRM	0.014	167.10>152.10	20.00	239.10>167.10	8.00	239.10>72.00	16.00
45-2 Fosfthiazole-2	MRM	0.014	195.00>109.00	10.00	195.00>60.00	22.00	195.00>139.00	6.00
45-3 Permethrin ethyl	MRM	0.014	304.10>168.10	12.00	318.10>165.10	12.00	318.10>182.10	12.00
45-4 Isopropalin	MRM	0.014	280.10>238.10	8.00	280.10>133.10	18.00	280.10>165.10	16.00
45-5 Isodrin	MRM	0.014	192.90>157.00	20.00	192.90>128.00	26.00	262.90>192.90	28.00
45-6 Cyprodinil	MRM	0.014	224.10>208.10	16.00	224.10>197.10	22.00	224.10>131.10	14.00
45-7 Isopropalin-methyl	MRM	0.010	199.00>121.00	14.00	241.10>121.10	22.00	0.00>0.00	0.00
45-8 MGK 264-2	MRM	0.014	164.10>98.00	12.00	111.10>82.00	8.00	164.10>67.00	8.00
45-9 (E)-Chlorfenvinphos	MRM	0.014	323.00>267.00	16.00	267.00>159.00	18.00	267.00>203.00	12.00
45-10 Dimethametryn	MRM	0.014	212.10>122.10	12.00	212.10>94.00	22.00	212.10>71.00	18.00

▮ LabSolutions Insight Heightening the Efficiency of Multi-Analyte Data Analysis

With LabSolutions Insight software, quantitative results for a complete series of data files can be displayed side-by-side for comparison and QC review. All of the chromatograms for a selected target compound can be displayed simultaneously, making it easy to review the detected peaks and confirm the quantitative results. Color-coded QA/QC flags quickly identify any outliers that require further examination.

■ Color-coded Quantitative Flags

In LabSolutions Insight, quantitative results can be compared to established criteria, and any outliers are color-coded for easy identification and further review. Five color-coded criteria levels can be defined, making it easy to determine which data points are outliers, and which specific QC criteria were not met. Any changes made to calibration curves or manual peak integration are immediately reflected in the color-coded flags.



■ Status Review Function


This function can be used to specify the status of all compounds and samples for their management. By specifying a status, the progress of data analysis work can be accurately recorded and reported.

#	Flags	Data Filename	Status	Sample Type	#	Name	Status
<input checked="" type="checkbox"/>		Conc-3_001	Pending	Standard	<input checked="" type="checkbox"/>	MPA-gluc	Pending
<input checked="" type="checkbox"/>		Conc-3_002	Accept	Standard	<input checked="" type="checkbox"/>	d3-MPA...	Pending
<input checked="" type="checkbox"/>		Conc-3_003	Review	Standard	<input checked="" type="checkbox"/>	MPA	Pending
<input checked="" type="checkbox"/>		Unknown-1	Rerun	Unknown	<input checked="" type="checkbox"/>	d3-MPA	Pending

Rerun
Pending Accept


■ Enhanced Accuracy Controls

Accuracy control results, such as retention time differences, calibration curve linearity and peak shapes (tailing), can be assessed visually. In addition, the QC chart function in LabSolutions Insight allows visual confirmation of variations in target compounds across multiple data sets. This is useful when evaluating variations in internal standard substances between samples, for example.



■ System Configurations Using Multiple Client Computers

Data acquired from multiple systems can be reviewed or confirmed using client computers connected via a LAN or other network. If multiple systems are used, data obtained from each system can be reviewed from any client computer. Even in the case of multiple analysts using the same system, the ability to separate analytical work from measurement work improves work efficiency.



File management on a file server is recommended for systems with more than five users.

GC-MS/MS Databases Allow Users to Start Analysis Immediately

Offering both quantitative and screening applications, these databases contain pre-optimized MRM transitions and collision energies required for GC-MS/MS analysis, allowing users to start an analysis immediately just by using the AART function to automatically modify the retention times.

Smart Database Series Supporting Accurate Quantitative Determinations in MRM Analysis

The series contains optimized analysis conditions, so users can start an analysis immediately without investigating the conditions. This database is for accurate quantitative determinations using standard samples.

Smart Pesticides Database

It covers the pesticides (530 compounds) subject to GC-MS analysis and used inside and outside Japan. The database also contains information on compounds that can be used as internal standards. Therefore, it also supports analysis with the internal standard method.



Smart Environmental Database

The database contains information on 527 compounds including polychlorinated biphenyl, brominated flame retardants, dioxins, polycyclic aromatic hydrocarbons, and organochlorine pesticides, as well as their stable isotope labeled compounds.



Smart Metabolites Database

The database contains 525 compounds including metabolites contained in blood, urine, cells and other biological samples. It also contains information on the stable isotopes of 22 major metabolites, which can be used as internal standards.



Smart Forensic Database

The database is registered with 201 forensic toxicological substances often involved in poisonings, such as drugs of abuse, psychotropic drugs, pharmaceuticals, and pesticides.



Quick-DB Series Enables Quick Screening Without Using Standards

Preregistered with calibration curve information from the internal standard methods, this database allows users to calculate quantitative values without analyzing standard samples. It can be used for screening when quick confirmation of quantitative results is required.


Quick-DB GC/MS Residual Pesticides Database

The database contains information on calibration curves created using pesticide surrogates as the internal standards, thus enabling the simultaneous screening of 491 residual pesticides without using standard samples. The registered calibration curves have been created by grouping target pesticides based on similar behavioral and physical properties, and allocating pesticide surrogates with similar physical properties to each group as internal standard substances. This allows highly accurate quantitative values to be calculated.



Allocation of Pesticide Surrogates Based on Physical Properties of Target Pesticides

Target pesticides




Pesticide surrogates

- Group 1
Dichlorvos-d
- Group 2
Acephate-d
- Group 4
Iprobenfos-d

Allocation

Sensitivity Variation Adjustments Using Pesticide Surrogates





Pesticide surrogates can also be used to adjust differences in target pesticide responses observed in calibration curves and a customer's instrument due to variations in GC-MS/(MS) sensitivity.

● Quick-DB Forensic Toxicology Database

This database contains information on sample preparation, data acquisition, and analysis of 68 compounds often involved in poisonings. Due to the simple QuEChERS method used for the sample preparation, even first-time users can pretreat samples easily. The performance control function, which automatically assesses the system status, quickly determines problematic areas. As a result, users can analyze toxicological substances while consistently maintaining the highest status level.

Using a Calibration Curve for a Pretreated Sample to Correct for Recovery Rates

By using the calibration curve information from a pretreated sample, accurate quantitative values can be obtained. The calibration curve information is used to correct differences in recovery rates for each target compound or correct for effects caused by adsorption, which occurs when creating calibration curves using standards.



If specific control parameter criteria are not satisfied, the optimal steps for restoring the system status are displayed in a message box.



Quick-DB Forensic Launcher Screen


● Off-Flavor Analyzer

This analysis system contains information on the major odor-causing substances identified from previous problems and associated sensory information, thus enabling the reliable identification of odor-causing substances. It allows users to calculate quantitative values easily without using standard samples, and to identify odor-causing substances through comparison with registered odor thresholds. The system also supports the use of a sniffer, enabling the efficient confirmation of odors using the predicted retention time display function. In addition, support is provided for a total system, including pretreatment units such as HS, SPME, and thermal desorption units using MonoTrap.

Compound Name (E)	Ref. Index 1	Comment (E)	Threshold
Benzophenone	2471	Almond, Burnt Sugar	10
2,4,6-Trichloroanisole	2801	Lardine	100
1-Tetralone	2182	Cinnamon	1000
gamma-Endosulfone	2784	Sweet, Pinen, Fruit	1
Octen-3-one	3002	Fruit	1

Primary odor components GC/MS analytical conditions Sensory information

Comparison of Concentration to Odor Threshold Values




Odor-causing substances can be identified through comparisons of odor thresholds with the concentrations of off-flavor components identified in a chromatogram.

Name	Conc	Unit	Threshold	Description
Benzophenone	2.543	pg/mg	10.000	Almond, Burnt sugar
2,4,6-Tribromophenol	2241.933	pg/mg	100.000	Lodoform

<p>Under Threshold Value</p> <p>Odor threshold value: 10.000 pg/mg</p> <p>Concentration: 2.543 pg/mg</p> <p>Benzophenone</p>	<p>Over Threshold Value</p> <p>Concentration: 2241.933 pg/mg</p> <p>Odor threshold value: 100.000 pg/mg</p> <p>2,4,6-Tribromophenol</p>
--	---









Since some odor components have a low odor threshold, low concentration levels need to be detected to identify the causative substances. High-sensitivity MRM/SIM analysis by GC-MS/MS can reliably detect even trace components near the odor threshold (a few pg/g).









Qualitative Analysis Using GC-MS/MS

GC-MS/MS allows users to freely select ions cleaved by electron ionization (EI), and then cleave precursor ions via collision induced dissociation (CID), thus enabling the detailed analysis of partial structures. A product ion scan can be used to easily discriminate structural isomers and regioisomers, which are hard to identify with an EI mass spectrum obtained via GC/MS. This mode is optimal for predicting partial structures.

Example of Structural Isomers




Example of Regioisomers




High-Speed Scanning Control (Advanced Scanning Speed Protocol, ASSP)

UFMS
ULTRA FAST MASS SPECTROMETRY



Newly Patented Technology (ASSP)



Black: 1,111 u/sec
Red: 5,000 u/sec
Blue: 10,000 u/sec

Equipped with a function to automatically optimize the rod bias voltage during high-speed data acquisition, the system can collect data stably through scans as well as product ion scans at 20,000 u/sec (patent: US6610979).

With the GCMS-TQ8040, ten or more product ion scans can be configured for a single measurement, and three or more user-specified analysis modes (Scan, MRM, product ion scan) can be combined.


	Compound Name	Start Time (min)	End Time (min)	Acq. Mode	Event Time(sec)	Scan Speed	Start m/z	End m/z	Precursor m/z	GE
1-1		2.30	43.00	Q3 Scan	0.050	10000	43.00	500.00		
1-2		2.30	43.00	MRM	0.050					
1-3	Amine m/z 244	2.30	43.00	Product Ion Scan	0.005	20000	20.00	45.00	44.00	15.00
1-4	Amine m/z 258	2.30	43.00	Product Ion Scan	0.005	20000	20.00	59.00	58.00	15.00
1-5	Amine m/z 272	2.30	43.00	Product Ion Scan	0.008	20000	20.00	73.00	72.00	15.00
1-6	Amine m/z 286	2.30	43.00	Product Ion Scan	0.008	20000	20.00	87.00	86.00	15.00
1-7	Amine m/z 298	2.30	43.00	Product Ion Scan	0.010	20000	20.00	99.00	98.00	15.00
1-8	Amine m/z 310	2.30	43.00	Product Ion Scan	0.010	20000	20.00	101.00	100.00	15.00
1-9	Amine m/z 312	2.30	43.00	Product Ion Scan	0.010	20000	20.00	113.00	112.00	15.00
1-10	Amine m/z 314	2.30	43.00	Product Ion Scan	0.010	20000	20.00	115.00	114.00	15.00
1-11	Amine m/z 326	2.30	43.00	Product Ion Scan	0.010	20000	20.00	127.00	126.00	15.00
1-12	Amine m/z 328	2.30	43.00	Product Ion Scan	0.010	20000	20.00	129.00	128.00	15.00
1-13	Amine m/z 340	2.30	43.00	Product Ion Scan	0.010	20000	20.00	141.00	140.00	15.00
1-14	Amine m/z 342	2.30	43.00	Product Ion Scan	0.015	20000	20.00	143.00	142.00	15.00

Inert Flow Path Achieves High-Sensitivity GC/MS Flow Path

Minimizing adsorption and other losses in the flow path from sample injection to the detector is important for stable, high-sensitivity measurements using GC/MS. The flow lines in the GCMS-QP series and the TQ series consist of high-quality, highly reliable consumable parts, so even trace concentrations of components can be detected with high sensitivity and favorable repeatability.

Glass liner

The Restek glass liner recommended for GCMS analysis uses a proprietary inactivation technology to dramatically suppress active sites. After packing into the insert, the wool is subjected to a complete inactivation treatment. This product is controlled throughout from production to final inspection to provide 100 % satisfaction.



Blue: Sky liner (with wool)
Red: A liner other than Sky liner (with wool)

Micro-syringe



Autosampler syringes feature improved durability, clarity, and accuracy, achieving reliable injection accuracy.

GC septum

Our lineup now includes low bleed septa, which maintain optimal seal performance even when the injection cycles are increased, and can be used even at high temperatures. This reduces sensitivity variations due to leaks.

Ferrules and gold gasket

The high-quality Vespel ferrule is easily attached and designed to resist leaking. The gold gasket is inactive, and adsorption does not occur.

Ion source




Designed with a shield that blocks radiant heat generated by the filament, and an ion source treated with an oxide coating, active spots inside the ion source are not prone to occur, which enables high-sensitivity analysis with long-term stability.

Capillary columns



For the SH-Rxi series, a high-quality fused silica like no other is used as the raw material. Our proprietary surface inactivation technology and optimal process to mask silanol groups result in a low-bleed column with very impressive inactivation performance, even with respect to polar compounds comparable to acidic and basic substances.



Optional Units Enable a Variety of System Configurations

▮ Various System Configurations


For GC-MS/MS analysis, different system configurations may be required depending on the application and sample-introduction needs. The GCMS-TQ8040 offers a wide variety of system configurations and sample-introduction devices to enable an expanded range of applications.



▮ Chemical Ionization and Negative Chemical Ionization

In addition to electron ionization (EI), chemical ionization (CI) and negative ion chemical ionization (NCI) are available for the GCMS-TQ8040. The CI method softly ionizes sample molecules and is effective for confirming molecular weights, while the NCI method enables selective, high-sensitivity detection of compounds with functional groups that have a large electron affinity such as halogens.

The Quick-CI function allows users to switch between the EI and CI modes with the same ion source, so EI and CI data can be obtained from the same sample. If confirmation of molecular ions is difficult with the EI mode, molecular information can be collected from the CI mode data.




▮ DI-2010 Direct Sample Inlet Device



Direct sample injection (DI) is a method in which a sample is injected directly into the ion source without passing through the gas chromatograph (GC). This is an effective method for measuring the mass spectra of synthetic compounds, and can be used easily for the usual GC-MS


configuration. By using this in combination with Quick-CI, EI and CI mass spectra can easily be collected.

Components that are thermally degradable or difficult to vaporize are not suited to GC analysis. Their mass spectra can be obtained easily using the DI probe.



The AOC-6000 supports three sample injection methods: liquid sample injection, headspace (HS) injection, and solid phase micro extraction (SPME) injection, so samples in a variety of forms can be analyzed. It allows the sample injection method to be switched automatically, enabling different sample injection methods to be combined in a continuous operation.

With the automatic syringe exchange and stirring function, standard samples can be prepared automatically with a variety of dilution levels, and everything from the creation of calibration curves to the quantitative determination of unknown samples can be fully automated.



The HS-20 headspace sampler provides strong backup for the analysis of volatile components at every stage from research to quality control departments. The high-sensitivity electronic cooling trap enables quantitative and qualitative determination of trace components that cannot be detected with a conventional headspace sampler.

The HS-20 transfer line is built into the GC unit, which makes it easy to combine the HS-20 with the AOC-20 liquid sample injector, as well as to switch between these units.




The OPTIC-4 multimode sample inlet is a GC injection port that enables a variety of sample injection modes for GC-MS, including large-quantity injection, inlet derivatization, thermal desorption, and DMI (difficult matrix introduction).

Combining this with an autosampler enables automatic replacement of inserts, improving productivity in multisample analyses.



DMI (Difficult Matrix Introduction) Mode Simplifies Pretreatment



An undiluted shampoo solution was placed directly in a microvial and measured in DMI mode, enabling the selective analysis of volatile components.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.